

## Introduces HT Universal Layering Ingots



Swiss NF Metals, Inc. is proud to introduce SNF PressCeram® (HT) High Temperature Universal Layering Ingots.

The SNF (HT) ingots allow the laboratories the freedom to press thin veneers in the range of .4mm to .8mm with the added benefit of using most High Fusing PFM ceramic materials, either to stain or layer, giving the ceramist the freedom to use existing materials and being cost effective.

SNF PressCeram HT Ingot is a high temperature stable feldspar material containing leucite Press-Glass ceramic with a CTE of 13.9.

The increased pressing temperature results in greater stability and a controlled coefficient of thermal expansion which allows trouble-free layering and characterization with the technicians chosen conventional High Fusing PFM material.

Indications for thin Veneers, \*Anterior All-ceramic crowns, Inlays and Onlays.  
*\*single units only*



MONARCH METALS  
For USA Call: 1.800.531.2444



1.800.554.5504